





01	Glass Liner Light-Curing Lining Material	09	Porcelain Edge 9.5% Hydrofluoric Acid Etching Ge
02	Edge Line Glass-Ionomer Lining Cement	10	Edge Prep 17% EDTA Cream
03	Edge Fill Glass-Ionomer Filling Cement	11	Edge Seal Resin Resin-Based Root Canal Sealer
04	Edge Com Nano-Hybrid Universal Composite	12	Edge Seal (MTA) Bioceramic Sealer (MTA-Based)
05	Edge Flow Nano-Hybrid Flowable Composite	13	Edge Seal (TCS) Bioceramic Sealer (TCS-Based)
06	Edge Bond Universal Bonding Materials 5 th Generation	14	SUPER-Files Endodontics Rotary NiTi Files
07	Edge Bond SE Self Etch Bonding Material 7th Generation	15	E-Flex Endodontics Rotary NiTi Files
08	Blue Edge 37% Phosphoric Acid Etching Gel	16	Edge File Endodontics Rotary NiTi Files

HAMERZ MEDICAL

Knowledge-Based Company



HAMERZ Medical Co., as a leading player in the dental healthcare industry, strives to become the first choice for dentists and specialists in this field.

Committed to excellence, HAMERZ Medical Co. leverages its strong infrastructure and a professional, experienced team to manufacture and deliver advanced and high-quality dental materials and equipment in the fields of endodontics, prosthetics, and restorative dentistry, adhering to the latest European and global industry standards and technology.

At HAMERZ Medical Co., our priority is the quality and safety of our products. We have established a robust and comprehensive quality management system, incorporating international standards such as ISO 13485, ISO 9001, and GMP. Our research and development team, comprised of expert professionals from reputable universities, focuses on innovation and product improvement, actively contributing to the development of new dental technologies.

We are committed to environmental preservation and offering eco-friendly products. In our production processes, we employ methods, materials, and equipment selected to reduce waste, enhance energy efficiency, and protect the environment.

Customers play a vital role in HAMERZ Medical Co.'s expansion of its product portfolio. We prioritize providing continuous, comprehensive services and support to our customers. Alongside our professional and responsive sales team, we aim to provide custom solutions and continuous enhancements to meet our customers' needs.

HAMERZ Medical Co. takes pride in covering a broad spectrum of the Iranian market, entering the global market, and reducing the final product's cost. Our mission is to facilitate dentists' access to the latest, advanced technology products. And our vision is to be recognized in this industry as a familiar and trustworthy name for dentists.

Glass Liner

Light-Curing Lining Material
Syringe **2.0 ml (3.6 g)**



Glass Liner is a one-component light-curing, flowable, and radiopaque glass-ionomer lining material. Glass Liner can be used as a lining under all types of filling materials, including composites, amalgams, indirect restorations, extended fissure sealing, and for the treatment of smaller lesions. Glass Liner is thixotropic and thereby easy to handle. It also contains fluorides which remineralize the tooth enamel.

Specification

Direct and easy injection
Effectively prevents secondary caries
High fluoride release and remineralization
High compressive strength and strong adhesion to dentine
Bacteriostatic effect and high biocompatibility
Sets well under all types of filling materials
Ideal for extended fissure sealing
Thixotropic properties
High radiopacity
Easy handling
Suitable for minimal lesions and cavity preparation

Technical Specification

Density: 1.7-1.8 g/cm³

Compressive Strength: 200 Mpa

Curing depth: 2 mm Polymerization time: 30 s Viscosity: 170-290 Pa.s

Shelf Life: 3 years (6 months after opening)

Edge Line is a two-component self-cure glass-ionomer cement with minimal shrinkage during hardening, results in less sensitivity compared to other composite materials. This product is suitable for lining under all kinds of filling materials, direct repair of small lesions and repairing deciduous teeth. It can be used as fissure sealant and shows excellent radiopacity. This product also contains fluoride preventing secondary caries.

Edge Line

Glass-Ionomer Lining Cement
Powder bottle **15.0g**Liquid Bottle **10.0ml**

Specification

Strong adhesion to both dentin and enamel
Remineralization with high fluoride release, preventing secondery caries
High compressive strength
Long-term biocompatibility
Transparent material without negative effect on the filling
Radiopaque
Self-cure

Technical Specification

Working Time: 100 s Setting Time: 180-220 s

Film Thickness: 10 micrometers Compressive Strength: 180 MPa

Water solubility: 0.2% Adhesion to dentin: 5 MPa

Shelf life: 3 years



Edge Fill

Glass-Ionomer Filling Cement

Powder bottle **15.0g**Liquid Bottle **10.0ml**

EDGE FILL Glass-lonomer Filling
Powder 159

Liquid 10m

Edge Fill is a two-component self-cure glass-ionomer cement for direct restorations of the cavity classes III and V, tunnel preparations, and deciduous teeth. It may also be used as fissure sealant. This product bonds well to the enamel and the dentine and can be polished after 15 min. Edge Fill is radiopaque and releases fluoride into the tooth. This product is available as a flexible powder/liquid system in two different shades (A2 and A3).

Specification

Rapid application compared to composite
Application in anterior and proximal contact areas
Radiopaque
Self-cure
Long-term biocompatibility
High compressive strength
Excellent adhesion to the dentin and the enamel
For recommended-shaped and abrasion defects, and in deciduous teeth

Technical Specification

Working Time: 90-120s Setting Time: 180 s

Film Thickness: 10 micrometers Compressive Strength: 180 MPa

Water solubility: 0.2% Adhesion to dentin: 4 MPa

Shelf life: 3 years

EDGE COM is a nano-hybrid light-cure composite for restorations of all cavity classes in the anterior and posterior areas. EDGE COM is easy to handle because it does not stick to the instruments and at the same time adapts prefectly to the cavity walls. The product shows a very low shrinkage. EDGE COM also generates esthetic results because it is high gloss polishable and translucent. Thanks to its high strength and outstanding esthetics, Edge COM is suitable for veneers.

Specification

Nano-hybrid Light cure

Repairing of all cavity classes and restorations in the posterior and anterior areas

Applicable for direct veneers

High strength and wear-resistance

The low shrinkage of only 1.7%, minimizing the risk of secondary decay

Thixotropic, not sticky to tools

Excellent polishability

Desired translucency

Radiopaque

Color stability

Technical Specification

Compressive strength: > 45 MPa Flexural strength: > 155 MPa Flexural modulus: > 6 GPa Shrinkage: 1.7% linear Curing Depth: 2 mm

Color Shades: in accordance with VITA® shade guide *

Ceramic filler: 77% wt Curing time: 20s

*VITA is the registered trademark of vita-zahnfabrik

Edge Com

Nano-Hybrid Universal Composite
Syringe **4.0 g**



www.hamerzmedical.com

Edge Flow

Nano-Hybrid Flowable Composite
Syringe **2.0 g**



www.hamerzmedical.com

EDGE FLOW is a flowable, light-curing composite for direct restorations of the cavity classes III, IV and V as well as repairs of bis-acrylate materials. EDGE FLOW is thixotropic, which makes it agreeable to handle. The composite does not run away from the place of application while is fluid when manipulated. EDGE FLOW is also nicely polishable and delivers esthetic restorations. The shrinkage is 1.7%, which reduces secondary caries.

Specification

Suitable for direct restoration of class III, IV, V caries
Filling small cavities of class I, II caries
Filling the anterior areas of the cavity-class III, IV
Bonding of ceramic and composite restorations
Extended fissure sealing for premolar and molar teeth
Low shrinkage of only 1.7 %, reducing the risk of secondary decay
Thixotropic
Light cure
Excellent polishability
Optimal viscosity
Radiopaque
Easy repair and very good handeling
Color stability

Technical Specification

Compressive strength: > 80 MPa Flexural strength: > 5.6 MPa Flexural modulus: > 80 MPa

Shrinkage: 1.7 % Curing depth: <2 mm

Color Shades: in accordance with VITA® shade guide *

Ceramic filler: 54% wt Curing time: 30-40 s EDGE BOND is a bonding agent for the total-etch technique with excellent adhesion to dentin and enamel. It ensures a secure adhesion between the tooth and composites. The material shows a good flowable consistency and is easy to use. EDGE BOND also adheres on glass-ionomer linings.

Edge Bond

Universal Bonding Materials 5th Generation Bottle **5.0 ml**

Specification

Excellent adhesion to dentin and enamel
Resistant to evaporation
Used in the connection between dentin and composite
High bonding strength (13.5 MPa)
Very high permeability
Water tolerant (insoluble in water)
No gaps in the margin areas
Adherance to lining materials including glass ionomer
Easy to use and long-term durability

Technical Specification

Tensile strength: 12 MPa

Notched-edge shear strength: 41 MPa

Curing time: 15-30 s Viscosity: 0.3-0.8 Pa.s



Edge Bond SE

Self Etch Bonding Material 7th Generation
Bottle **5.0 ml**

EDGE BOND Self-Etch (SE) is a self-etching, light-curing one-component bonding for adhesive restorations. The material contains methacrylate. One step application of EDGE BOND SE makes it time saving material.



One step etch and bond without prior etching
Elastic connection between tooth and composite
Perfect bonding of composites to dentin, enamel and lining materials, including glass ionomer
Very high permeability
Contains methacrylate
No gaps in the margin areas
Hydrophilic
Easy to use and long-term durability
Water tolerant (insoluble in water)



Bonding strength: 13.5 MPa Viscosity: >0.5 dynamic Pa.s

Curing duration: 10 seconds with 800mW/cm² power



Blue Edge

37% Phosphoric Acid Etching GelSyringe **3.0 ml**

BLUE ETCH is in the form of 37% phosphoric acid gel. Blue Etch creates microscopic spaces in the tooth enamel (increasing the surface roughness) which results in very good bonding (micro-mechanical retention). When drying, the etched surface appears white. Blue Etch also removes

When drying, the etched surface appears white. Blue Etch also removes surface contamination and helps keep tooth enamel moist. (It removes the smear layer and the hybrid layer.)

Specification

37% phosphoric acid solution Gluten and glycerin free Total-Etch system function Increasing the bonding strength

Suitable for use before composite, bonding, adhesive and all kinds of sealants Ideal viscosity and the ability to permeation of all gaps

Not spreading vertically and to unwanted parts

Unique self-limiting chemical characteristic: 2 μm Average depth of penetration within 15 to 20 seconds

Removing quickly and easily with water spray

The possibility of precise deployment in areas with severe access restrictions

Technical Specification

Etching time: 30 sec for enamel and 15 sec for dentin Depth of etching: 2 μm Thixotropic



Porcelain Edge

9.5% Hydrofluoric Acid Etching GelSyringe **3.0 ml**



PORCELAIN EDGE (containing 9.5 % hydrofluoric acid) is a hydrofluoric acid, which is specifically used to create an etched surface for broken porcelain restorations and all-porcelain restorations. It can also be used to etch ceramic veneers before connecting orthodontic brackets.

Specification

Hydrofluoric acid gel 9.5 %

Can be used for: surface etching of veneers, inlays, porcelain and ceramic onlays, broken teeth inside the mouth

Preparation of porcelain surface for stronger bonding with resin and composite Can be washed quickly and easily with water spray

The possibility of precise deployment in areas with severe access restrictions

Technical Specification

Etching time: 60 sec Depth of etching: 2 µm

Thixotropic

Edge Prep

17% EDTA Cream Syringe **5.0 ml**

EDTA 17% Solution is buffered to a neutral pH and is an effective calcium binding or chelating agent used to facilitate instrumentation of root canals and for smear layer removal. It decalcifies the canal walls making it easier to enlarge and shape the canal with files and reamers.

Specification

EDTA solution is used to wash the tooth canal
A clear, colorless and odorless liquid
Softens the tooth root and dissolves calcium with the chelate effect
Easy to open and clean the calcified canal
Cleans the dentin wall after removing the smear layer
Removes inorganic substances
Neutral pH
In the form of cream, gel or liquid

Technical Specification

Peroxide containing



Edge Seal Resin

Resin-Based Root Canal Sealer
Syringe 10.0 g



EDGE SEAL RESIN is a resin based root canal sealer material with excellent sealing ability and excellent biocompatibility. This product is a non-staining material to the resin. Edge Seal shows a good radiopacity and easy handling. As a pre-mixed product, the application is very fast and easy.

Specification

More effective obturation in shorter time
Provides a complete, thorough seal
Great radiopacity
Provides syringe delivery to the apical third
Low shrinkage and no gap between the sealer and the canal wall
Low solubility and resistant to breakage over time
Available as double syringes (base and catalyst) with a ratio of one to one
Suitable for side canal filling
Great biocompatibility
Creating a thin layer around gutta percha
Ease of automatic mixing of base and catalyst
No discoloration on teeth
Excellent flowability
Insoluble in tissue fluids
Suitable for sensitive pulps

Technical Specification

Flowability: 22 mm Setting time: 50 h (at 37 °C) Film Thickness: 20 µm

Solubility: 0.9%wt

EDGE SEAL (MTA) is a paste-type root canal sealer that has excellent physical and biological properties of MTA. It is premixed and pre-loaded in a syringe that allows direct application of the sealer into the root canal without requiring powder/liquid mixing. The product has outstanding flowability, which makes it possible to completely fill the root canal system including accessory and lateral canals. It is eugenol-free and will not impede adhesion inside the root canal.

Edge Seal (MTA) is developed for root canal filling and root perforation repair, effectively preventing periapical tissue irritation as well as secondary infection.

Specification

High biocompatibility
Low shrinkage, superior dimensional stability
Hydrophilic
Bonding the hydroxyapatite with the dentin wall
High antibacterial effect
Homogeneous particle size distribution
Fast setting time
Easy and fast root canal filling
Increased formation of hard tissue
Excellent film thickness
Ready-to-inject paste
Complete apical sealing

Technical Specification

Fast setting time: 10-12 minutes
Easy retrieval with NiTi files
Hermetic apical sealing
High radiopacity
Promotion of hard tissue formation
No postoperative pain after overfilling

Edge Seal (MTA)

Bioceramic Sealer (MTA-Based)
Syringe **3.0 g**



www.hamerzmedical.com

Edge Seal (TCS)

Bioceramic Sealer (TCS-Based)

Syringe 3.0 g



EDGE SEAL (TCS) have high tissue compatibility. This product has a very low shrinkage showing superior dimensional stability. Edge Seal (TCS) is hydrophilic and leads to bonding the hydroxyapatite with the dentin wall.

Specification

Containing tricalcium silicate
High biocompatibility
Low shrinkage, superior dimensional stability
Hydrophilic
Bonding the hydroxyapatite with the dentin wall
High antibacterial effect
Homogeneous particle size distribution
Fast setting time
Easy and fast root canal filling
Increased formation of hard tissue
Excellent film thickness
Ready-to-inject paste
Complete apical sealing
Good regenerative ability

Technical Specification

Fast setting time: 10-12 minutes
Easy retrieval with NiTi files
Hermetic apical sealing
High radiopacity
Promotion of hard tissue formation
No postoperative pain after overfilling

The SUPER FILE is an economic endodontics rotary instrument made of high quality controlled memory Nickel Titanium wire with cutting edges and passive tip. SUPER FILE, having progressive taper comes in two colors, blue and gold, provides good flexibility and fracture toughness used to mechanically shape and prepare the root canals during endodontics therapy or to remove the root canal filling material while performing re-treatment. Super-File also shows good resistance against torsional fracture and cyclic fatigue failure, making it suitable for anterior and posterior teeth. Having safe tips, SUPER FILE avoids canal transportation and ledge formation.

SUPER-Files

Endodontics Rotary NiTi Files

Specification

Gold system

Treatments of anterior teeth with moderate to mild curved root canals Treatments of posterior teeth with partially easy root canals Higher cutting efficiency compared to the blue one

Blue system

Treatment of complex root canals such as curved and narrow and calcified ones Higher flexibility compared to the gold one

Technical Specification

HV-200 micro-Vickers hardness: >365 Torsional fracture strength: 10.9 mN.m

Fracture degree: >490 Stiffness: 8.0 mN.m



E-Flex Files

Endodontics Rotary NiTi Files

E-Flex is a highly flexible endodontics rotary instrument specifically designed to negotiate curved canals shows excellent resistance against torsional fracture and cyclic fatigue failure. E-Flex in two colors, gold and blue, comes in different tapering systems, 0.04, 0.06 and progressive for varying canal anatomies and a great all-round file system. E-Flex provides excellent cutting efficiency and high flexibility, makes it a good choice to reduce ledge formation, apical extrusion and canal transportation. Using different techniques of manufacturing and having variable pitch, E-Flex prevents locking and brings safe treatments.

Specification

Gold system

Treatments of anterior teeth with moderate to mild to moderate curved root canals Treatments of posterior teeth with partially easy root canals Higher cutting efficiency compared to the blue one

Blue system

Treatment of complex root canals such as curved and narrow and calcified ones Higher flexibility compared to the gold one



Technical Specification

HV-200 micro-Vickers hardness: >268 Torsional fracture strength: 10.9 mN.m

Fracture degree: >490 Stiffness: 8.0 mN.m

EDGE FILE is an outstandingly flexible endodontics rotary instrument using very high quality NiTi alloy which improves strength and flexibility and makes it suitable for difficult canals. Proprietary heat treatment processing makes the file extremely flexible which reduces the bounce back effect compared to the other NiTi rotary files. In addition, our special electro-polish process increases sharpness and strength of the file which leads to high efficiency in cutting and removing debris. Employment of different tapering systems results in prevention of locking and bringing high quality safe treatment and re-treatment.

Specification

Optimum flute diameter for minimally invasive preparation Maximizes file cutting efficiency due to a very specific cross section and non cutting tip design Increasing sharpness and strength of file using special electro-polishing technique Reduced handle length for increased posterior access

Technical Specification

HV-200 micro-Vickers hardness: >370 Torsional fracture strength: 10.9 mN.m

Fracture degree: >490 Stiffness: 8.0 mN.m

Edge Files

Endodontics Rotary NiTi Files



www.hamerzmedical.com



















Block 36, Piece 1/2, 4th Street, Payam Economic Free Zone Alborz, IRAN. | Tel: +98 26 34008673 | www.hamerzmedical.com Postal Code: 3187411512





